



*Patient Warming
System* **PW-11**

New generation of body temperature management

Hypothermia, when the core body temperature becomes lower than 36 °C can cause complications such as increased oxygen consumption, adverse cardiovascular events, coagulation dysfunction, decreased efficacy of anesthetic drugs, delayed anesthesia recovery, and delayed postoperative recovery.

Acer Biomedicals patient warming system can accurately and efficiently manage the patient's body temperature, to help prevent perioperative hypothermia, reduce intraoperative volume of blood needed in cases of bleeding, reduce postoperative complications, and promote faster postoperative recovery.

Patient Warming System Features

FastPre-heat- Only 20 seconds from turn on to ready for use

Accurate Control- temperature range from 33 °C -44 °C with adjustable increments of 0.5 °C
 High Quality blanket design- Soft comfortable blanket with even air-flow distribution. Tear-proof and Water-Proof

Safe and reliable- disposable blankets for infection prevention.

Product Comparison

	Patient Warming System	Competitors
Convectionairheating	+	—
LCDDisplay	+	—
AdjustableTemperatureRange (0.5°Cincrements)	+	—
AdjustableAir-Flow(5Levels)	+	—
ImportedFan(Germany)	+	—
Noiselevel<50dB	+	—



ApplicationRange

Commonly used in anesthesiology, ICU, emergency room, disaster etc.

RegisteredAddress: Ground Floor Plot no-183 Pocket N, Sector-2, DSIIDC Industrial Area, North West, Delhi, 110039

WorkAddress: Industrial Unit-39, Block-B, Sector-63, Noida, Gautam Buddha Nagar, Uttar Pradesh-201301

Web: www.accuremedical.in, **Contact:** +91 8077818101, 0120 4916914

Email: info@accuremedical.in, ampl2207@gmail.com

Tollfree No.: 1800-419-9029