



ACCUDEP-22

Standard Configuration:

Manual defibrillation, AED, 3/5-lead ECG, RESP, Thermal Recorder

Optional:

Pacer, NIBP, PR, EtCO₂, SpO₂

Safety Standards:

Physical Characteristics

Size:	295mm×252mm×316mm
Weight	5.2kg (Including 1 battery)
Screen Size:	7/8.4" TFT screen
Resolution	800 × 480
Waveforms:	Max 4 waveforms

Operation Environment

Temperature:	0~45
Humidity:	°C 10%~95%, non-condensation
Atmosphere Pressure:	700hPa~1060hPa
Ingress Protection:	IP44
Power requirement:	100-240V~,50/60Hz±3Hz
Battery type:	Rechargeable Lithium-ion battery
Battery capacity:	7500mAh, d.c.14.8V 5000mAh, d.c.14.8V

Battery number:	1
Battery recharging Time:	7500mAh Battery: Less than 2 hours to 80% and less than 3 hours to 100% with equipment power off 5000mAh Battery: Less than 1.5 hours to 80% and less than 2.5 hours to 100% with equipment power off

Battery backup:	Monitoring Mode: no less than 6 hours Defib Mode: 210 times (360J charge at intervals of 1minute without recording);
-----------------	---

Brightness:

Indicator

Pacing Mode: 4.5 hours (Load:50Ω, frequency: 80bpm, current: 60mA, without recording)
5000mAh Battery:
Monitoring Mode: no less than 4 hours
Defib Mode: 120times (360J charge at intervals of 1minute without recording);
Pacing Mode: 3hours (Load:50Ω, frequency: 80bpm, current: 60mA, without recording)
Manual from X to 100, X refers to the darkest brightness (X is 10 by default)

Two alarm indicators
Power indicator
Battery indicator
Maintain indicator
QRS beep and alarm sound
Operating key sound

Interfacing

USB interface
RJ45 interface
AC power input
Multi-functional connector

Date storage

Alarm Event:	200 groups
Patient profiles:	100 groups
Patient Events:	1000 groups
Wave Review:	10 min
NIBP Review:	2000 groups
Trend Graph:	160 hours
Trend Table:	160 hours
Voice recording:	Max 240 min in total;

Marked events	(Up to 60 min for each patient)
Power-off storage:	Available
Alarm:	Yes
	User-adjustable High and Low 3-level Limits;
	Prioritized audible and visual alarm
Network:	Connected to Central Monitoring System by hardwire/wireless

Recorder

Type:	Built-in; Thermal array
Channel:	Max 3 channel waveforms
Real-time recording:	3s, 5s, 8s, 16s, 32s, Continual
Speed:	6.25mm/s, 12.5mm/s, 25mm/s, 50mm/s
Record width:	50mm
Resolution:	8dot/mm (Horizontal and vertical)

Background grid:	Configurable
External printer:	Yes

Defibrillation

Operating mode:	Manual Mode, AED Mode, , Synchronous Defibrillation
Waveform:	Biphasic truncated exponential waveform, with impedance compensation
Defibrillation pathway:	External defibrillation
Electrode type:	External defibrillation paddles, multifunctional electrode
External defibrillation electrode paddles:	Supports charging, discharging and energy selection; Charging completion indicator
Charge Time: (Battery power)	Less than 3 seconds to 200 Joules with a new, fully charged battery Less than 7 seconds to 360 Joules with a new, fully charged battery
Charge Time: (AC power)	Less than 4 seconds to 200 Joules; Less than 8 seconds to 360 Joules
Energy accuracy:	±1.5J or ±10% of setting, whichever is greater, while 50Ω impedance ±2J or 15% of setting, whichever is greater, while 25Ω, 75Ω, 100Ω, 125Ω, 150Ω, 175Ω impedance
Patient Impedance Range:	20~300Ω(External defibrillation);

Defibrillation proof:	Type CF: ECG, RESP, SpO2, NIBP, PR;
	Type BF: EtCO2

Manual Mode

External defibrillators:	1J~360J, 25 types (1/2/3/4/5/6/7/8/9/10/15/20/30/50/70/100/120/150/170/200/220/250/270/300/360J)
Synchronous Cardioversion:	Energy transfer begins within 60ms of the R wave from internal Sync signal Energy transfer begins within 25ms of the External Sync signal

AED

Output Energy:	Adjustable: 100-360J
Number of electric shocks	Adjustable: once, twice, 3 times
Types can be AED:	VF & VT
AED maximum time required for cardiac rhythm analysis to be ready for discharge:	Battery power supply: 18s AC power supply: 21s

Noninvasive Pacing

Waveform:	Monophasic square wave pulse
Pulse Width:	20ms or 40ms
Accuracy:	±5%
Pacing Mode:	On-demand or fixed
Pacing frequency:	30 ppm to 210 ppm
Accuracy:	±1ppm or ±1.5% (whichever is greater)
Pacing output:	0 mA to 200 mA
Accuracy:	±5% or ±5mA, whichever is greater
Speed-down pacing:	Pacing pulse frequency reduced to 25% of original value.

Monitoring

ECG (leads)

Lead Type:	3 leads ECG, 5 leads ECG, AUTO
Lead selection:	5-lead: I; II; III; aVR; aVL; aVF; V 3-lead: I; II; III

Multi-lead synchronization analysis:	Available
--------------------------------------	-----------

ECG sensitivity:	Auto, 1.25 mm/mV (×0.125), 2.5 mm/mV (×0.25), 5 mm/mV (×0.5), 10 mm/mV (×1), 20 mm/mV (×2), 40 mm/mV (×4),
------------------	--

Accuracy:
Sweep speed:

Accuracy: Less than $\pm 5\%$
Heart Rate: 6.25 mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s
Less than $\pm 10\%$

Alarm limit range

Adult: 15~300bpm
Pediatric: 15~350bpm
Accuracy: $\pm 1\text{bpm}$ or $\pm 1\%$ (whichever is greater)
Adult:
High limit: (low limit+2bpm) ~ 300bpm
Low limit: 15bpm~ (high limit-2bpm)
Pediatric:
High limit: (low limit+2bpm) ~ 350bpm
Low limit: 15bpm~ (high limit-2bpm)

Resolution:

Accuracy: 1 bpm

Bandwidth: $\pm 1\text{bpm}$

Monitoring: 0.5~40Hz (-3.0dB~+0.4dB)
Diagnosis: 0.05~150Hz (-3.0dB~+0.4dB)
Surgery: 1~20Hz (-3.0dB~+0.4dB)
ST: 0.05~40Hz (-3.0dB~+0.4dB)

CMRR:

Monitoring: $> 105\text{dB}$
Diagnosis: $> 90\text{dB}$
Surgery: $> 105\text{dB}$
ST: $> 105\text{dB}$

Input Impedance:

$\geq 5\text{M}\Omega$

Input signal range:

$\pm 8\text{mV}$

HR trigger threshold

200 μV

Lead off detection

Measuring electrode: $< 0.1\mu\text{V}$

current:

Driving electrode: $< 1\mu\text{V}$

Pacemaker pulse

Manual selection when the

suppression switch:

pacemaker is turned on

Analog output:

Magnification: 1:1000;
Accuracy: $\pm 5\%$
Bandwidth: 0.5Hz ~ 40Hz
Delay: $\leq 35\text{ms}$

ST Detection:

-2.0mV~+2.0mV (-20.0mm~+20.0mm)
0.01mV
-0.8mV ~ +0.8mV: $\pm 0.02\text{mV}$ or $\pm 10\%$;

Resolution:

Accuracy:

ST analysis review

System noise:

Calibration voltage

Arrhythmia Analysis:

Pacemaker detection:

ECG (paddle)

Lead Type: Heart Rate measurement & alarm range:

Resolution:

Accuracy:

Others: Unspecified

20 groups

Less than 25 μV

1 mV; Accuracy: $\pm 5\%$

26 Types

Detectable

Single lead ECG

Adult: 15~300bpm

Pediatric: 15~350bpm

1 bpm

$\pm 1\%$ or $\pm 1\text{bpm}$ (whichever is greater)

Defib: 1~20Hz (-3dB~+0.4dB)

Defib: $> 105\text{dB}$

Bandwidth:

CMRR:

Input Impedance:

$\geq 5\text{M}\Omega$

Input signal range:

$\pm 8\text{mV}$

HR trigger value

200 μV

Arrhythmia Analysis:

5 Types, ASY, VF, VT, PNC, and PNP

Respiration

Method:

Thoracic Impedance Method

RR measurement

Adult: 0~120bpm

range:

Pediatric: 0 ~ 150bpm

Accuracy:

7~150bpm: $\pm 2\text{bpm}$ or $\pm 2\%$ (whichever is greater)

0~6bpm: unspecified

Adult: 10s~60s Ped: 10s~40s

Apnea Alarm:

$\pm 5\text{s}$

Accuracy:

Audible and visual alarm; alarm

Alarm:

events reviewable

NIBP

Method

Automatic oscillometric

Work mode:

Manual / Automatic/Continuous

Interval Time:

Adjustable

1/2/2.5/3/4/5/10/15/30/60/90/120/180/240/480/720 min

Continuous: 5min

Adu/Ped: 120s

Maximum

measurement cycle

Measurement Unit:

mmHg / kPa selectable

Pressure types

Systolic, Diastolic, Mean

Range of systolic

Adult Mode: 5.3~36kPa

pressure

(40~270mmHg)

Pediatric Mode: 5.3~26.7kPa

(40~200mmHg)

Range of diastolic pressure:	Adult Mode: 1.3~28.7kPa (10~215mmHg) Pediatric Mode: 1.3~20kPa (10~150mmHg)
Range of mean pressure:	Adult Mode: 2.7~31.3kPa (20~235mmHg) Pediatric Mode: 2.7~22kPa (20~165mmHg)
Over pressure protection:	Adult: 39.6kPa (297mmHg) Pediatric: 32kPa (240mmHg) Tolerance: ± 0.4 kPa (± 3 mmHg), if ± 0.667 kPa (± 5 mmHg), if exceeds the above range, the monitor can still display normally, but the accuracy is not considered
Accuracy:	
Alarm limit:	Same as the range of measurement
PR from NIBP:	40~240bpm
Resolution:	1bpm
Accuracy:	$\pm 3\%$ or ± 3 bpm, whichever is greater

SunTech NIBP

Regulatory compliance:	YY 0670-2008
Initial inflation range:	Adult: 16~37.3kPa (120~280mmHg) Pediatric: 10.7~22.7kPa (80~170mmHg)
Maximum measurement cycle:	Adult: 130s Pediatric: 90s
Over pressure protection:	Adult/Pediatric: 40.0kPa (300mmHg)
Static pressure measurement range:	0kPa~40.0kPa (0mmHg~300mmHg)
Resolution:	
Range of systolic pressure:	

Range of diastolic pressure:	± 0.4 kPa (± 3 mmHg) Adult: 5.3~34.7kPa (40~260mmHg) Pediatric: 5.3~21.3kPa (40~160mmHg)
Range of mean pressure:	Adult: 2.7~26.7kPa (20~200mmHg) Pediatric: 2.7~16kPa (20~120mmHg)
PR from NIBP	Adult: 3.5~29.3kPa (26~220mmHg) Pediatric: 3.5~17.7kPa (26~133mmHg) 30~220bpm

Accuracy: $\pm 2\%$ or ± 3 bpm, whichever is greater

Nellcor SpO2

Measurement range: 0~100%
Resolution: 1%
Accuracy: $\pm 2\%$ (70~100%, Adu/Ped, non-motion)
1~69% unspecified
Alarm range: 20~100%
PR Measurement
Range: 20~300bpm
Resolution: 1bpm
Accuracy: ± 3 bpm (20~250bpm)
Unspecified (251~300bpm)
Alarm range: 20~350bpm

MASIMO SpO2

Measurement & alarm range 1~100%
Resolution: 1%
Accuracy: $\pm 2\%$ (70~100%, Ped/Adu, non-motion)
 $\pm 3\%$ (70~100%, motion);
1~69% unspecified
Alarm range 1~100%
PR Measurement
Range 25~240bpm
Resolution: 1bpm
Accuracy: ± 3 bpm (non-motion)
 ± 5 bpm (motion);
Alarm range: 20~350bpm
PI value: 0.02~20%
Resolution: 0.01% (0.02~9.99%)
0.1% (10~20%)
SIQ: Available

SpO2

0~100%
Measurement & alarm range:

Resolution:

Accuracy: 1%
 $\pm 2\%$ (70~100%, Ped/Adu, non-motion)
0~69% unspecified

PR Measurement

Range: 20~254bpm
Resolution: 1bpm
Accuracy: ± 2 bpm

Alarm range: 20~350bpm
 PI value: 0.05~20%
 Resolution: 0.01% (0.05%~9.99%)
 0.1% (10.0%~20.0%)

Accuracy:
 SIQ: unspecified
 Available

MASIMO EtCO₂ (Sidestream)

Measurement range: 0~190mmHg, 0~25vol%
 (at 760mmHg)
 Accuracy: Standard environment 22±5°C,
 1013±40kPa:
 a) 0~114mmHg:
 $\pm(1.52\text{mmHg} + \text{reading} \times 2\%)$
 b) 114~190mmHg: not defined
 All environment:
 a) 0~114mmHg:
 $\pm(2.25\text{mmHg} + \text{reading} \times 4\%)$
 b) 114~190mmHg: not defined
 1mmHg or 0.1% or 0.1kPa
 Resolution: 0~150rpm
 awRR range: ±1rpm
 awRR accuracy:

Response time: <3 s

Respironics EtCO₂ (Sidestream)

Measurement range: Loflow:
 0~150mmHg, 0~19.7%, (0~20kPa)
 (at 760mmHg)
 CapnoTrak:
 0~99mmHg, 0~13.03%, 0~13.2kPa
 (at 760mmHg)
 Loflow:
 Accuracy: ± 2mmHg (0~40mmHg)

± 5% of reading (41 – 70mmHg) ±
 8% of reading (71 – 100mmHg)
 ±10% of reading (101~150mmHg)
 (In 25°C, if RR>80rpm, accuracy is
 12% of reading

CapnoTrak:
 ±2mmHg (0~38mmHg)
 ±10% of reading (38~99mmHg)
 RR influence to EtCO₂
 (0~99mmHg):
 -2~0.5mmHg (0-40bpm)
 (-6% of reading)~0.5mmHg (41-
 70bpm)
 (-14% of reading)~0.5mmHg
 (71~100bpm)

Resolution: 1mmHg
 Loflow: 2~150rpm
 CapnoTrak: 0, 2~100rpm
 awRR range: ±1rpm
 awRR accuracy: